

# Patient chair Malin 4



## GENERAL

Malin 4 patient chair is popular in several treatment areas because it combines good ergonomics and a high functionality. The stable design with large wheels and robust mechanical design provides safety for both patient and staff.

Malin 4 is electrical and easy to adjust in height with the hand control. Malin 4 is available in different models and with plenty of accessories to customize the chair. The chair has a manual brake bar on the back wheels, electric brake is available as an option.

### MALIN 4 ROCK

120-00103-00



Malin 4 Rock is our appreciated base model that can be outfitted with accessories for different examination situations. Seat and back rest can be tilted to get the desired angle.

Lifting capacity	150kg / 331lb
Weight	41kg / 90lb
Height adjustment	540-840mm / 21-33in
Wheel frame	530x580mm / 21x23in
Backrest angle	-8° to +20°
Seat angle	+4° to -10°
Battery operated	

### MALIN 4 SWING

120-00103-01



Malin 4 Swing has a swivel seat that can be rotated  $\pm 90^\circ$ , to enable the patient to comfortably be seated from the side with the chair parked. The seat is conveniently turned to the correct position by the device table while the base of the chair remains in position.

Lifting capacity	150kg / 331lb
Weight	41kg / 90lb
Height adjustment	540-840mm / 21-33in
Wheel frame	530x580mm / 21x23in
Battery operated	

### MALIN 4 SLIM

120-0103-02



Malin 4 Slim has a narrower wheelbase in the front which makes it easier to manoeuvre the chair in places with limited space, e.g. close to device tables with only one lifting column.

Lifting capacity	150kg / 331lb
Weight	41kg / 90lb
Height adjustment	540-840mm / 21-33in
Wheel frame	530x580mm / 21x23in
Backrest angle	-8° to +20°
Seat angle	+4° to -10°
Battery operated	

## STANDARD SEAT COLOURS



### **140-00135-00 blue vinyl**

Vinyl has an airy, ventilated slick texture.

Seat 450x470mm / 18x19in, Backrest 360x410mm / 14x16in



### **140-00247-02 black art leather**

Artificial leather is more durable.

## ACCESSORIES

There are a variety of accessories that all have been developed in close cooperation with medical professionals, to optimize the chair for different treatments and to get an ergonomic and comfortable environment for both patients and staff.

### **140-00128-01 Extended lifting height, 100mm/ 4in extra**

### **140-00126-02 Wireless hand control**

### **140-00126-03 Wireless hand control for electric brake**

The IR hand control comes with a strap and bracket/holder to be placed where it's suitable; on the chair or by the working station.

### **140-00250-00 Electric brake**

One button for both brake and release on either the fixed or IR hand control



### **140-00315-00 Wheels with directional control**

Makes the chair easier to move around when back wheels have a directional lock



### **150-00102-00 Armrest standard left**

### **150-00102-01 Armrest standard right**

### **150-00103-00 Armrest collapsible left**

### **150-00103-01 Armrest collapsible right**



### **140-00129-03 Foot support, rubber (not for Malin slim)**

The foot support is raised along with the chair and can be adjusted in height or be folded up.




### **150-00132-00 Driving handle 450mm / 18in**

### **150-00132-02 Driving handle 600mm / 24in**



### **190-00018-00 Head support**





rini

Rini Ergoteknik AB  
Fågelsångsvägen 3B  
SE-186 42 Vallentuna  
Sweden  
Phone +46 8 594 77170  
Email [info@rini.se](mailto:info@rini.se)  
[www.rinicompany.com](http://www.rinicompany.com)



rini  
MedTech

Rini Medtech Inc  
USA  
Phone +1-331-230-6711  
Email [info@rinimedtech.com](mailto:info@rinimedtech.com)  
[www.rinimedtech.com](http://www.rinimedtech.com)

